Inkontinenztherapie-Stimulator

BioBravo

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Content

Please read the Operators Manual carefully before using the device and pay attention to important references !

| | Page |
|----------------------------------------------------------|-------------|
| I. Incontinence Therapy BioBravo | |
| Delivery & Accessories - Measurement & Weight - EC | 23 |
| Introduction - What is BioBravo ? - Indication | 24 |
| Contraindication | 25 |
| Advice on Application of Electrodes | 26 |
| Technical Data - Safety Adjustments | 27 |
| Maintenance - Warranty | 28 |
| II. Display Elements of the BioBravo | |
| Display and Operating Elements | 29 |
| Description of the Device | 30 |
| III. Application of the BioBravo | |
| Settings of The Preset Programmes | 31 |
| Settings of the Individual Programmes | 32-33 |
| IV. Therapy for Incontinence | |
| What is Incontinence ? | 34 |
| Different Types of Incontinence | 35 |
| Programmes 01 - 09 | 36-40 |
| V. Control of Therapy | 41 |
| VI. Literature | 42 |
| VII. Waste Disposal | 43 |
| VIII. Copyright | 43 |

Delivery & Accessories

Included in the delivery:

010E-003: 1 device BioBravo, 2 lead wires, 4 self adhesive electrodes (40 x 40 mm), 4 self adhesive electrodes (40 x 80 mm), two 9 V batteries and 1 Operators Manual

010E-038: 1 device BioBravo, 2 lead wires, 1 vaginal probe (IncoProbe V), 1 lubricant (75 g), two 9 V batteries and 1 Operators Manual

010E-039: 1 device BioBravo, 2 lead wires, 1 rectal probe (IncoProbe R), 1 lubricant (75 g), two 9 V batteries and 1 Operators Manual

Accessories:

011E-006: 1 unit self adhesive electrodes 40x40 (4 pieces)

011E-007: 1 unit self adhesive electrodes 40x80 (4 pieces)

011E-051: 1 unit self adhesive electrodes 32mm Ø (4 pieces)

011E-008: 1 vaginal probe “IncoProbe V”

011E-009: 1 rectal probe “IncoProbe R”

011F-012: 2 lead wires

011F-013: 2 pieces 9V batteries

011F-037: 1 tube electrode jelly (75g) for probes

011E-047: 1 operators manual BioBravo

011E-050: 1 storage box

BioBravo may only be used with original accessories.

Measurement & Weight

Measures of the device: H 14,1 cm / W 7 cm / D 4 cm

Weight of the device: 180 g (incl. battery)

Electromagnetic Compatibility

The device conforms to the electromagnetic compatibility requirements of European Standard EN 60601-1-2:2001 “Medical Electrical Equipment, Part 1: General Requirements for Safety - Collateral Standard: Electromagnetic Compatibility - Requirements and Tests”

Introduction – What is BioBravo ?

BioBravo is a neuromuscular stimulator, which has been developed according to latest findings by MTR+ in consultation with the medical profession and patients. This device is a modern 2-channel stimulator, which is often used with vaginal- or rectal-probes. Additionally this unit allows the use of 2 pairs of skin-electrodes (= 4 electrodes) to work simultaneously, too. **BioBravo** is a stimulation device, which allows both, high-efficient muscle stimulation for the pelvic muscles as well as being suitable for the calming of the hyper nervous bladder. This device can therefore be used for all forms of urinary and/or stool incontinence without any foreign assistance. Another program tries to make the activation of the peristalsis possible, so that this device can be used for remedying of obstruction, too.

Additionally the three freely configurable programmes enable the specialist to compile specially tailored programmes for his own individual purpose. This characteristic of the **BioBravo** offers therefore each physician or physiotherapist, but also the experienced patient, the possibility of a focused coordinated treatment for the respective state of health. As these programmes, after storing, can be changed at any time, for example to be altered according to the progress of therapy, an almost unrestricted variation range is available with this device to the user.

After termination of the therapy, **RehaBravo** can be re-used by other patients. Please clean it before with a mild disinfectant cleaner.

Indication

BioBravo was designed to aim at the following:

Therapy of incontinence:

- Elimination of light to mid strength faecal incontinence (PRG 7, 8)
- Elimination of light to mid strength urinal incontinence
 - Stress-incontinence (PRG 3, 4)
 - Urge-incontinence (PRG 1, 2)
 - Mixed-incontinence (PRG 5, 6)
- Activation of peristalsis (PRG 9)

Contraindication

BioBravo may under no circumstances be used in the following cases:

- Patients with a pace maker as functional disturbances can occur
- Patients with myocard damages (damage to the cardiac muscle) or arrhythmia (inconsistent heart beat)
- When operating machinery with a potential of danger
- Patients with epilepsy
- When inflammation of the vagina, anus or urinary tract is present
- Patients with raised urine remainder
- Patients, who have previously suffered from urine retention

Please pay unconditional attention to the following ground rules:

- Do not stimulate directly through ostesynthesis material (metal in the body)
- Do not stimulate transcerebral (via the head)
- Do not stimulate via the Carotis-Sinus-Nerve (carotid artery)
- Do not stimulate in the area of the Sinus knot (heart)
- Do not stimulate on the Bulbous (eye)
- Do not stimulate in the Larynx-Pharynx-Area (larynx /pharynx)
- Do not stimulate in places with skin defects
- Do not position the electrodes directly on the spine

Take special care with the use of the electric stimulation:

- If after an operation the healing process can be interrupted through a muscle contraction.
- After acute bleeding or bone fractures with a tendency to bleed.
- In the case of the menstruation.
- If in the case of damage to the perception nerves the normal sensation in the skin has been lost.
- If at the same time the connection is made with a surgical high frequency apparatus, this could result in burn marks underneath the electrodes.
- If device operation occurs close to short waves or a microwave, discrepancies of the initial values of the stimulator could occur.
- The application of the electrodes nearby the thorax can increase the risk of ventricle cilia.
- In the case of pregnancy, stimulation should only be applied after consultation with the general practitioner.
- If in the case of treatment with the incontinence probe, pain or irritation appear in the area of the pelvis, or bleeding in-between menstruation, interrupt the therapy and consult your doctor.
- Should skin irritation occur whilst using the **BioBravo**, please interrupt the application and inform your doctor, so that an allergic reaction can be ruled out.
- Patients with implants should stimulate only after a previous consultation of an experienced doctor.

Please, keep this device stored out of the reach of children.

Advice on Application of Electrodes

Please empty the urinary bladder and clean the area to be treated.

Make sure that BioBravo is turned off (no display on the LCD display).

Now connect the single poled ends of the lead wires with the vaginal or rectal probe, respective the skin-electrodes.

Then fix the output plug(s) of the lead wires into channel A or B of BioBravo.

Furthermore make sure that all plugs fit completely into the electrode sockets as well as the device sockets.

If you are using a probe, make it gliding capable with water or a gel suitable for electrotherapy.

Please don't use any oil, ointment or cream, since all these have an isolating effect.

Now place the probe or the electrodes according to the instructions of your doctor or physiotherapists.

If you should use a probe, please clean it after use carefully with a normal detergent and dry it completely.

Electrodes/probes which were damaged due to an improper handling may not be used any more.

Electrodes/probes with a damaged cable isolation (e.g. visible copper wire may not be used any more).

Electrodes/probes may not be changed (e.g. by means of scissors or other tools).

Do not pull the electrodes/probes by force on the electrode cable.

Apply the electrodes on clean and unwounded skin only.

In the case of strong bodily hair growth, the hair should be shortened with scissors.

Neither stretch the skin nor the electrodes when applying the electrodes.

To remove the electrodes from the skin, please lift the electrodes at the rim towards the middle.

Effective current densities of more than 2 mA/cm² require an increased attention of the user with regard to the intensity adjustment. All customary electrodes which comply with the requirements set out in the Council Directive 93/42/EEC relating to the Medical Device Directive can be used, if the connections match the cables of **BioBravo**. We recommend the use of our vaginal- (IncoProbe V) or rectal- (IncoProbe R) probes.

Technical Data

| | |
|--------------------|-----------------------------------------------------------|
| Channels: | Two isolated independent channels |
| Waveform: | Asymmetrical, rectangular, bi-phasic with zero DC current |
| Intensity: | 0 - 90 mA |
| Frequency: | 2 - 200 Hz |
| Width of impulse: | 50µs - 450µs |
| Ramp up/down time: | 0,1 seconds - 9,9 seconds |
| Time setting: | variable: 1 minute - 9 hours + 59 minutes |
| Electrical supply: | 9 volt block battery |
| Resistor: | 500 Ohm 1000 Ohm 2200 Ohm |
| U max: | 48 V 80 V 120 V |

Environmental conditions for storage, transportation & use: -10 to +50 degree centigrade, 0-90 % atmospheric humidity, air pressure of 700-1060 hPa

Explanation of Symbols


Attention!

Please consider
Operators Manual

0123
CE-Sign

Grade of protection BF

This device is quality controlled and fulfils the rules and regulations set by the German law for medical products (MPG) as well as the guidelines of the European Community 93/42/EEG for medical products. Therefore this device carries the CE-sign „CE 0123“.

You can find the serial number of this device in the battery compartment.

Manufacturer: MTR+ Vertriebs GmbH, Kamenzer Damm 78, 12249 Berlin

Safety Adjustments

1. The +buttons get inactive after one minute running time of a program to avoid an unintentional increase of intensity . This lock is switched off again, only by pressing the +button three-times repeatedly. Now an increasing of intensity is possible again for another minute.
2. If the electrical circuit should be interrupted for example by losing an electrode or by a broken lead wire, the device switches off immediately, if it has been above an intensity of 6 mA.
3. If any phase comes to its end the intensity falls down to 8 mA automatically and must be regulated by the user on a pleasant measure again.

Maintenance

BioBravo has been designed to be maintenance free. Nevertheless, please adhere to the following advice:

If needed please clean **BioBravo** and the probes with a mild detergent, water and a damp cloth.

Do not expose **BioBravo** to extreme humidity and do not submerge under water.

You may clean the cables with a damp cloth. You may also from time to time treat the cables with talcum powder in order to prevent cracks and to further ongoing usage.

If over a longer period **BioBravo** is not in use, please take the battery out of the device.

Please subject the device to technical safety control in regular intervals in accordance with the local regulations.

Warranty

MTR+ Vertriebs GmbH provides a 24 month warranty, commencing from the date of our delivery (invoice date), that in case of failure, the **BioBravo** will be serviced or replaced free of charge. All returns must first be authorised by MTR+ Vertriebs GmbH in advance.

Excluded from this warranty are failures of the **BioBravo** which result from excessive use, careless or improper handling as well as damage by force.

The warranty loses its lawful meaning if the producer MTR+ Vertriebs GmbH has not given its consent in writing prior to repairworks or servicing being carried out. Without this written consent the manufacturer will not return the device.

Display and Operating Elements

„ON/OFF“: By pressing the button once, the device will either be switched on or off.

„+/- buttons“: This button regulates the intensity for the channels A and B. The device also has two independently controllable output channels A and B. When pressing the respective arrow buttons the intensity will either be increased or decreased. 60 seconds after the last intensity regulation, a further increase may not be carried out due to an electronic barrier. But if the „+“ button is pressed three times, the barrier is released and the intensity can be altered again. A reduction of the intensity is possible at any time.



Whilst configuring the programmes 10 – 12 one switches between the programmable parameters by pressing the left „+/-“ button (channel A), and through pressing the right „+/-“ button (channel B) the desired value/data for respective parameter can be set. In this way, the application time, the stimulation, the pause, rising pulse and falling pulse, as well as synchronous or asynchronous stimulation can be selected.

„PRG“: With this button the preset programmes 01 – 09 or the three free programmable programmes 10 – 12 can be selected for individual use. In addition, this button stores the current adjustment in the programming mode and ends it.

„SET“: With the SET-button one can select the different programme phases. With every pressing of this button, the phase in progress is terminated the next programme phase is started. This applies to either the running of the programmes 01 – 09 and also to the programming and the sequence of programmes 10 – 12. If one keeps the SET-button pressed for 3 seconds whilst being in programmes 10 – 12, the programming part of the programme is then entered.

Description of the Device

Settings of the Preset Programmes 01–09

Preparations:

Make sure that **BioBravo** is switched off. Cleanse the area of the skin which is to be treated with clear water. Following that, make sure that all plugs are completely connected to the respective electrode sockets as well as device sockets. Place the vaginal-/rectal-probe or the adhesive electrodes onto the area that has been prior cleansed. The clip on the back of the unit allows to fasten the device at the clothes.

BioBravo is operated by an 9 volt block battery. Open the battery compartment at the back by pushing down the lid and put in the battery by hand. Please ensure the correct polarity in accordance with the markings. **BioBravo** is now operational. If the battery voltage is too low, the battery-symbol in the LCD-display will flash. Please replace batteries immediately. Empty batteries constitute special waste and therefore have to be disposed of at an appropriate collecting point. Never dispose of batteries via domestic or household waste.

Adjustment:

1. Press once the button „**On/Off**“.

BioBravo is now switched on. The program used last, appears on the LCD-display.

2. By pressing the programme button „**PRG**“ you can choose from the programmes 01 – 09 (pre-programmed).

3. This is followed by pressing the buttons „+/-“.

The selected programme now runs automatically. If the button „+“ is pressed repeatedly the intensity of stimulation increases. If the button „-“, is pressed repeatedly the intensity of stimulation is reduced. For the setting of the intensity it is vital that one feels the electric current but that pain is never incurred. **Please take into account that for safety reasons the intensity falls back at every changing of phase on the value of 8 mA and must be regulated by you on the desired standard again.** Each phase end is indicated by 5 peep tones.

4. If you wish to interrupt the set programme, then press the button „**PRG**“ once. The programme will run again when pressing the button „+“. In this case, the intensity must be adjusted newly, too.

5. **BioBravo** is switched off when pressing the button „**On/Off**“ or automatically one minute after the program is finished.

Settings of the Individual Programmes 10-12

1. Switch on **BioBravo** with “**On/off**” button.
2. Choose the desired program (13-15) with “**PRG**” button.
3. Press “**SET**” button for at least three seconds. **BioBravo** changes to the programming mode and the left watch, which indicates phase duration, flashes to the sign that this function can be changed.

In the following please use the left “+” button always then, if you want to store your choice and to change to the next step. The right “+/-” buttons are respectively meant for the change of the chosen and flashing parameter.

4. Select with the right “+/-”buttons the desired time for the first phase.
5. By pressing the left “+”button you will be able to store your choice and change to the next step. In the display now the mode flashes, which you adjust according to your wishes by pushing the right “+/-”buttons.
If you have chosen “**CONT**” or “**BURST**” you can determine frequency and after storage and change with the left +button respectively the pulse width by pushing the right “+/-”buttons.
If you have chosen “**WORK/REST**”- rhythm you can determine the following parameters one after the other by pushing the right “+/-”buttons.
 - frequency (**Hz**)
 - pulse width (**µ S**)
 - working time (**sec**)
 - rest period (**sec** = recovery time)
 - increase time of the impulse (**ramp up**)
 - way down time of the impulse (**ramp down**)
 - synchronous (**SYN**) or alternating (**ALT**) work of the channels A and B. If you have selected “**SYN**” you can program **delay**-time between the channels A and B.

If you have chosen “**MOD**” you can change the following parameters one after the other.

- low limit of frequency (**Hz down**)
- upper limit of frequency (**Hz up**)
- low limit of pulse width (**μ S down**)
- upper limit of pulse width (**μ S up**)
- low limit of intensity (the display right below shows values from 0.5 to 1.0). This means that the intensity of current can be selected between 50% (=0.5) of the original adjusted mA-value 100% (=1.0). The upper limit of the intensity, that is 100%, is always the mA value adjusted by the user originally. By the fact that the intensity, even in a modulated program, can never increase above the 100%-value it is ensured, that there will be not more current, then the user has individually regulated himself during the whole stimulation process.
- modulation time (possible choice: 5 sec. up to 60 sec, except intensity: pre-set 1 sec.)

6. By pushing the “**SET**”button, you reach the second phase of the program. You can analogously program these to the first sequence now. This way you can make up to 5 phases of your own in each of the three configurable programs. If you should want to reduce a configured program of 5 phases to less phases, you simply put the time to zero in the unwanted program sections.

7. The programming is completed by pressing the “**PRG**” button and the program shall be stored as long as you don’t want to change programming.

Abbreviations:

BURST: The burst mode delivers 2 blocks of 9 impulses each per second

CONT: In the continuous mode a permanent unchanged current takes place.

MOD: In the modulated mode the frequency and/or the pulse width and/or the intensity changes in a predefined rhythm automatically.

SEC: Second

Hz: Hertz is the measurement unit for the frequency of impulses per second.

μ S: Microsecond is the measurement unit for the length of every single impulse.

MA: Milli-amp is the measurement unit for the electrical current intensity.

SYN: Synchronous function of the channels A and B

ALT: Alternating function of the channels A and B

DELAY: Delay between channels A and B

Therapy for Incontinence

What is Incontinence ?

One uses the term incontinence in the case of spontaneous and uncontrolled loss of urine or stool.

**Incontinence can be treated !
Anyone stricken with incontinence can be helped !**

Incontinence is a widely spread problem. It is assumed to be realistic when stating that in many countries half of the people ageing 50 or more are affected by it. Even younger people can be afflicted with it.

Many folds are the causes for incontinence. Numerous factors - often combined – are responsible for it, such as malfunctioning of hormones, poor muscular functioning, injuries etc.

As a result of these varying causes, it is vital for an afflicted person to consult a doctor. The doctor alone will be able to initiate a suitable therapy.

The electrical stimulation is not stressful, is simple in its application and allows wide ranging use for treating bladder and discharge malfunction as well as faecal incontinence.

Forms of incontinence, upon which medicines have little effect are well therapeutible by electrical stimulation very often.

The aim of the electrical stimulation is as follows:

- focused training of the pelvic muscles
- strengthening of the bladder closing muscle
- attenuates the highly active bladder muscle
- regaining the physiological sensation of the bladder
- strengthening of the sphincter muscle
- activation of peristalsis

Different Types of Incontinence

BioBravo can be applied with the following types of incontinence.

Stress incontinence means that the group of muscle that belongs to the closing mechanism of the bladder exit is not fully functionable. When pressure rises (e.g. through sneezing, coughing, laughing or physical exertion) the loss of urine will occur even without micturition. The cause is a weakening of the pelvic muscle which can result from either weakness of the connective tissue, birth or prostate surgery.

Urge incontinence means that the closing mechanism of the bladder and urethra are mostly functioning. Often it comes to extrude movements of the bladder, which cannot be suppressed through mind power. One speaks of a hyperactive bladder muscle. Micturition is uncontrollable and the pressure of the bladder increases to the point, that the resistance of the healthy muscular closing system will be overcome. Possible causes are disorder of the superior nervous system.

A **mixed form of stress and urge incontinence** means that the cause lies with the weakness of the muscular closing system as well as a hyperactive bladder muscle.

Faecal incontinence means a disorder of the function of the sphincter anus externus. A controlled discharge of bowel contents is not possible anymore. The causes for the functioning disorder of the sphincter anus externus are manifold.

Programmes

| Progr.: 01 | | Urge incontinence 1 | | | | |
|-----------------|---------|---------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 20 | | | | |
| Mode | | Cont | | | | |
| Frequency work | Hz | 10 | | | | |
| Frequency rest | Hz | | | | | |
| Pulse duration | μ S | 250 | | | | |
| Modulation time | secs | | | | | |
| Ramp up time | secs | | | | | |
| Ramp down time | secs | | | | | |
| Work time | secs | | | | | |
| Rest time | secs | | | | | |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 20 min | | | | | |

Calming of the hyperactive bladder.

| Progr.: 02 | | Urge incontinence 2 | | | | |
|-----------------|---------|---------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 20 | | | | |
| Mode | | W/R | | | | |
| Frequency work | Hz | 10 | | | | |
| Frequency rest | Hz | | | | | |
| Pulse duration | μ S | 350 | | | | |
| Modulation time | secs | | | | | |
| Ramp up time | secs | 1 | | | | |
| Ramp down time | secs | 1 | | | | |
| Work time | secs | 6 | | | | |
| Rest time | secs | 15 | | | | |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 20 min | | | | | |

Calming of the hyperactive bladder.

Programmes

| Progr.: 03 | | Stress incontinence 1 | | | | |
|-----------------|--------|-----------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 20 | | | | |
| Mode | | W/R | | | | |
| Frequency work | Hz | 30 | | | | |
| Frequency rest | Hz | | | | | |
| Pulse duration | µS | 200 | | | | |
| Modulation time | secs | | | | | |
| Ramp up time | secs | 0,8 | | | | |
| Ramp down time | secs | 0,8 | | | | |
| Work time | secs | 5 | | | | |
| Rest time | secs | 8 | | | | |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 20 min | | | | | |

Training the weakened pelvic muscle. This program can be used with faecal-incontinence as well.

| Progr.: 04 | | Stress incontinence 2 | | | | |
|-----------------|--------|-----------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 10 | 10 | | | |
| Mode | | W/R | W/R | | | |
| Frequency work | Hz | 30 | 40 | | | |
| Frequency rest | Hz | | | | | |
| Pulse duration | µS | 450 | 300 | | | |
| Modulation time | secs | | | | | |
| Ramp up time | secs | 1 | 1 | | | |
| Ramp down time | secs | 1 | 1 | | | |
| Work time | secs | 6 | 6 | | | |
| Rest time | secs | 15 | 15 | | | |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 20 min | | | | | |

Training the weakened pelvic muscle. This program can be used with faecal-incontinence as well.

Programmes

| Progr.: 05 | | Mixed incontinence 1 | | | | |
|-----------------|--------|----------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 3 | 10 | 5 | 4 | 3 |
| Mode | | W/R | W/R | W/R | W/R | W/R |
| Frequency work | Hz | 3 | 10 | 20 | 30 | 40 |
| Frequency rest | Hz | | | | | |
| Pulse duration | µS | 250 | 250 | 250 | 200 | 200 |
| Modulation time | secs | | | | | |
| Ramp up time | secs | 0,8 | 0,8 | 0,8 | 0,7 | 0,7 |
| Ramp down time | secs | 0,8 | 0,8 | 0,8 | 0,7 | 0,7 |
| Work time | secs | 4 | 4 | 4 | 4 | 4 |
| Rest time | secs | 4 | 4 | 4 | 6 | 6 |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 25 min | | | | | |

Treatment of mixed urethral incontinence. Both the overactive bladder and the weakened pelvic muscle are treated.

| Progr.: 06 | | Mixed incontinence 2 | | | | |
|-----------------|--------|----------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 10 | 10 | 10 | | |
| Mode | | W/R | Cont | W/R | | |
| Frequency work | Hz | 12 | 10 | 35 | | |
| Frequency rest | Hz | | | | | |
| Pulse duration | µS | 400 | 200 | 300 | | |
| Modulation time | secs | | | | | |
| Ramp up time | secs | 1 | | 1,5 | | |
| Ramp down time | secs | 1 | | 1 | | |
| Work time | secs | 6 | | 5 | | |
| Rest time | secs | 15 | | 12 | | |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 30 min | | | | | |

Treatment of mixed urethral incontinence. Both the overactive bladder and the weakened pelvic muscle are treated.

Programmes

| Progr.: 07 | | Faecal incontinence 1 | | | | |
|-----------------|--------|-----------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 20 | | | | |
| Mode | | W/R | | | | |
| Frequency work | Hz | 40 | | | | |
| Frequency rest | Hz | | | | | |
| Pulse duration | µS | 200 | | | | |
| Modulation time | secs | | | | | |
| Ramp up time | secs | 1 | | | | |
| Ramp down time | secs | 1 | | | | |
| Work time | secs | 6 | | | | |
| Rest time | secs | 15 | | | | |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 20 min | | | | | |

Training the sphincter muscle. This program can be used with stress-incontinence as well.

| Progr.: 08 | | Faecal incontinence 2 | | | | |
|-----------------|--------|-----------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 5 | 10 | 5 | | |
| Mode | | W/R | W/R | W/R | | |
| Frequency work | Hz | 40 | 35 | 30 | | |
| Frequency rest | Hz | | | | | |
| Pulse duration | µS | 300 | 400 | 450 | | |
| Modulation time | secs | | | | | |
| Ramp up time | secs | 1 | 1 | 1 | | |
| Ramp down time | secs | 1 | 1 | 1 | | |
| Work time | secs | 5 | 6 | 5 | | |
| Rest time | secs | 12 | 12 | 12 | | |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 20 min | | | | | |

Training the sphincter muscle. This program can be used with stress-incontinence as well.

Programmes

| Progr.: 09 | | Activation of peristalsis | | | | |
|-----------------|---------|---------------------------|---------|---------|---------|---------|
| | | Phase 1 | Phase 2 | Phase 3 | Phase 4 | Phase 5 |
| Phase time | min | 10 | 10 | | | |
| Mode | | W/R | W/R | | | |
| Frequency work | Hz | 8 | 12 | | | |
| Frequency rest | Hz | | | | | |
| Pulse duration | μ S | 450 | 400 | | | |
| Modulation time | secs | | | | | |
| Ramp up time | secs | 1 | 1 | | | |
| Ramp down time | secs | 1 | 1 | | | |
| Work time | secs | 5 | 6 | | | |
| Rest time | secs | 10 | 12 | | | |
| Alternating | | | | | | |
| Synchronous | | | | | | |
| Overall time | 20 | | | | | |

In this program 09 it is tried to activate the peristalsis. To reach this target one uses both channels of the BioBravo in connection with self-adhesive electrodes, which are fixed crosswise on the body in the intestines area. As soon as the impulse appears, that the bowels should be emptied, the program can be stopped immediately.

- BURST = burst-mode
CONT = continuous mode
MOD = modulating mode (frequency and/or pulse range change automatically in accordance with a pre-given rhythm.
W/R = intermittent operating and resting mode

Modulation time = exponential

For example: Modulating time 10 seconds at a pulse range of 250 – 150 μ s means:

Start at 250 μ s – exponential decrease (first fast then slow) down to 150 μ s in five seconds, then increase (first slow then fast) back to 250 μ s. This complete cycle lasts for 10 seconds.

Control of Therapy

BioBravo has a memory store, where the information for time in use and the average mA current used for both channels can be read. The device summarizes the total time of use, up to a maximum of 999 hours and 59 minutes.

This way it is possible to see whether the intensity used by the patient has been high enough to get a treatment success, or if maybe intensity has been too strong, so that it was able to produce unpleasant side effects.

„LOCK“: If you want to lock the device, please open the battery compartment on the back of the unit. Press the concealed “LOCK”-button using the end of the lead wire until you hear a double beep. On the LCD-display appears “L:” and the symbol of the watch. By pressing the “PRG”-button you set the time to zero and then the device adds up all following application times. If you should have the desire to close all other programs for the use additionally, please press the “+/-”-buttons after the double beep. This way you add “P” to the display, so that you can see “L: P + watch-symbol”. After the storage of this choice with the “PRG”-button the device is locked now and adjustments cannot be changed anymore. To unlock the unit press the concealed switch once again and you will hear a single beep. The information for time in use and the average mA current used for both channels can be read. By pressing the “PRG”-button you bring the device to the Home position and the stored values are deleted.

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In accordance with European regulations [1], used electrical and electronic equipment may not be given to the unsorted waste anymore. The symbol of the litter bin on wheels indicates the necessity of separate collection.

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Reserved for changes.

Printing or copy of this manual - also in an abridged version - is forbidden and can only be allowed with a written approval of “MTR+ Vertriebs GmbH”.

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