

Incontinence-therapy



... the + of certainty



- 9 pre-set programs
- 3 additional five-phase custom programs
- Two channels
- Time setting from 1 min. to 10 hours
- Time count down monitored on LCD screen
- Automatic switch off after 4 minutes if unit is not in use
- Quick view to the most important treatment data or detailed control of each session of the last 60 days of treatment
- Real Time Clock (RTC)
- High safety level to protect the user
- CE 0123

BioBravo can be used for the following types of incontinence:

- **Stress Incontinence:** The group of muscles in the urethra and bladder which prevent urine from leaking through stress related situations are weakened. Sneezing, coughing, laughing and physical exertion cause urine to be lost because the muscles begin to stop functioning effectively. BioBravo stimulation helps recover and improve muscle response and strength.
- **Urge Incontinence:** The closing mechanism of the bladder and urethra function well. The bladder feels the urge to release urine with certain movements, positions or most of the time. The bladder muscle becomes hyper active due to the urge to urinate and weakens to the point where it cannot hold urine and leakage occurs.
- **Mixed forms of Stress and Urge incontinence:** The urethra and bladder do not function well and leakage occurs due to stress and a hyper-active bladder muscle.
- **Faecal Incontinence:** A disorder which effects the function of the sphinkter anus externus. A controlled discharge of bowel contents becomes very difficult or impossible.

Technical Data:

- Dual channel: individually isolated circuits
- Waveform: Asymmetrical, rectangular, bi-phasic with zero DC current
- Amplitude: 0 – 90 mA , Type: Constant Current
- Pulse width: 50 μ S – 450 μ S
- Pulse rate: 1 – 200 Hz
- Ramp up/down time: 0.1 – 9.9 seconds

1 Set BioBravo contains:

- 1 BioBravo
- 2 lead wires to connect 2 pairs of electrodes
- 4 batteries type AA
- 1 operators manual
- 1 carrying case and
- 4 self-adhesive skin-electrodes (40 x 40 mm)
- 4 self-adhesive skin-electrodes (40 x 80 mm)
= REF-Nr. 010E-003 or
- 1 vaginal probe IncoProbe V + 1 electrode-jelly
= REF-Nr. 010E-038 or
- 1 rectal probe IncoProbe R + 1 electrode-jelly
= REF-Nr. 010E-039