

## Multifunctional-stimulator



- 12 programs including 1 fully customizable 5 phase program
- Dual channel
- Timer and milliamps displayed on the LCD screen
- Custom program treatment times from 1 minute to 10 hours
- Quick view to the most important treatment data or detailed control of each session of the last 60 days
- Real Time Clock (RTC)
- High safety level to protect the user
- CE 0123

### Myolito combines three therapies in one device:

For the treatment of pain, four different operating modes can be used:

- Continuous flow of current
- Modulated Pulses
- Burst mode
- Han mode

For the rehabilitative training of muscles there are programs included that allow either targeted muscle growth and invigoration or that may serve for relaxation and detoxification. Muscle imbalances can be prevented and regeneration processes are accelerated.

The incontinence therapy includes treatments for both stress and urge incontinence as well as for faecal incontinence. Mixed incontinence can be treated easily with individually designed different phases in the custom program.

### Technical Data:

- Dual channel: individually isolated circuits
- Waveform: Asymmetrical, rectangular biphasic with zero DC current
- Amplitude: 0 – 80 mA - Type: Constant Current
- Pulse width: 50 µS – 450 µS
- Pulse Rate: 1 – 150 Hz
- Ramp up time: 0.1 - 9.9 seconds

### 1 Set Myolito:

- 1 Myolito Device
  - 2 Lead Wires to connect 2 pairs of electrodes/probe
  - 1 Battery 9V
  - 1 Operators Manual
- and
- 4 self-adhesive skin-electrodes (40 x 40 mm)
  - 4 self-adhesive skin-electrodes (40 x 80 mm)
  - = REF-Nr. 010E-205
- or
- 1 vaginal probe IncoProbe V + 1 electrode-jelly
  - = REF-Nr. 010E-206
- or
- 1 rectal probe IncoProbe R + 1 electrode-jelly
  - = REF-Nr. 010E-207