

Muscle Stimulation (EMS)



... the + of mobility

- 15 pre-set programs, incl. 1 program in “Dense-Disperse-Mode” according to Prof. Han
- 3 additional five-phase custom programs
- Two channels
- Time setting from 1 min. to 10 hours
- Time count down monitored on LCD screen
- Automatic switch off after 4 minutes if unit is not in use
- Quick view to the most important treatment data or detailed control of each session of the last 60 days of treatment
- Real Time Clock (RTC)
- High safety level to protect the user
- **CE 0123**

RehaBravo is designed to:

- Increase and strengthening of weakened and /or healthy muscles
- Increase in muscle endurance
- Prevention or delay of muscle degradation (atrophy) as a result of injury
- Maintain and increase range of motion
- Improved blood flow to the muscles
- Faster muscle recovery after exercise

Technical Data:

- Dual channel: individually isolated circuits.
- Waveform: Asymmetrical, rectangular, bi-phasic with zero DC current.
- Amplitude: 0 – 90 mA, Type: Constant Current.
- Pulse width: 50 μ S – 450 μ S
- Pulse rate: 1 – 200 Hz
- Ramp up/down time 0.1 – 9.9 seconds

1 Set RehaBravo: (REF-Nr. 010E-002)

- 1 RehaBravo
- 2 lead wires to connect 2 pairs of electrodes
- 4 self-adhesive electrodes (40 x 80 mm)
- 4 batteries type AA
- 1 operators manual
- 1 carrying case