

## Multifunctional Stimulator



... the + of quality in life



- 31 pre-set programs for:
  - muscular stimulation (EMS) – pain therapy (TENS)
  - incontinence therapy – sports/fitness/regeneration
- 3 additional five-phase custom programs
- Digital two channel-device with 4 exits for simultaneous use of up to 4 pairs of electrodes
- Alternating or synchronous work of the channels
- Time setting from 1 min. to 10 hours
- Time count down monitored on LCD screen
- Quick view to the most important training data or detailed control of each session of the last 60 days.
- Real Time Clock (RTC)
- High safety level to protect the user
- CE 0123

### MyoBravo combines several functions in one system:

The focus is on muscle stimulation. Here are programs included, which are used for specific muscular rehabilitation, as well as those that can help to improve performance of athletes or that can accelerate regeneration in sport.

Pain management programs are available in three different modes (continuous, modulated, burst), which can be usefully complemented by a fourth program, which contains a combination of two different modes for both goals of TENS therapy - blockade of nerves and increase of endorphins – in an automatic sequence.

The incontinence treatment includes multiphase programs for stress- and urge-incontinence as well as for faecal- and mixed-incontinence.

The custom programs allow to realize the own ideas regarding all these areas very comfortable. Of course, any custom-configured program can be changed anytime. The process of programming is very easy even without any experience.

### Technical Data:

- Dual channel: individually isolated circuits
- Waveform: symmetrical, rectangular, bi-phasic with zero DC current
- Amplitude: 0 – 90 mA , Type: Constant Current
- Pulse width: 50  $\mu$ S – 300  $\mu$ S
- Pulse rate: 2 – 200 Hz
- Ramp up/down time: 0.1 - 9.9 seconds

### 1 Set MyoBravo:

- 1 MyoBravo
- 4 lead wires to connect 4 pairs of electrodes/probe
- 4 batteries type AA
- 1 operators manual
- 1 carrying case and
- 4 self-adhesive skin-electrodes (40 x 40 mm)
- 4 self-adhesive skin-electrodes (40 x 80 mm)  
= REF-Nr. 010E-004 or
- 1 vaginal probe IncoProbe V + 1 electrode-jelly  
= REF-Nr. 010E-014 or
- 1 rectal probe IncoProbe R + 1 electrode-jelly  
= REF-Nr. 010E-015