

## Muscle-stimulator (EMS)



- 5 programs including 1 fully customizable 5 phase program
- Dual channel
- Timer and milliamps displayed on the LCD screen
- Custom program treatment times from 1 minute to 10 hours
- Quick view to the most important treatment data or detailed control of each session of the last 60 days
- Real Time Clock (RTC)
- High safety level to protect the user
- CE 0123

### Rehalito is designed to:

- Increase and strengthening of weakened and/or healthy muscles
- Increase in muscle endurance
- Prevention or delay of muscle degradation (atrophy) as a result of injury
- Maintain and increase range of motion
- Improved blood flow to the muscles
- Faster muscle recovery after exercise

### Technical Data:

- Dual channel: individually isolated circuits.
- Waveform: Asymmetrical, rectangular bi-phasic with zero DC current
- Amplitude: 0 – 80 mA - Type: Constant Current
- Pulse width: 50  $\mu$ S – 450  $\mu$ S
- Pulse Rate: 1 – 150 Hz
- Ramp up time: 0.1 - 9.9 seconds

### 1 Set Rehalito: (REF-Nr. 010E-201)

- 1 Rehalito Device
- 2 Lead Wires to connect 2 pairs of electrodes
- 4 self-adhesive electrodes (40 x 80 mm)
- 1 Battery 9V
- 1 Operators Manual